

Standard Autoinjector **Manufacturing System**

With ATS SuperTrak CONVEYANCE™ and Illuminate™ Software



Key Benefits

- Pre-engineered for lower design costs and faster delivery
- Suitable for any common three-piece autoinjector (additional part assembly options available)
- · Capable of handling external drive units
- · Accommodates both round and cut flange syringes (full vision inspection option)
- Various component loading options available from manual to fully automated
- Optimized compact floor space design
- Standardized documentation package
- · Low maintenance and minimum power/air usage

Key Facts

- Scalable Output Rate: 20-120 ppm
- MOEE ≥ 95%
- Reject Rate ≤ 1%
- Cleanroom Class.: ISO 8
- Footprint: 8m x 5m (for 120 ppm)
- Manufacturing Intelligence: Illuminate Base Module (OEE **Monitoring**)

Contact us for more information about our standard and custom solution offerings

ATS, SUPERTRAK CONVEYANCE, and ILLUMINATE are trademarks of ATS Automation Tooling Systems Inc. and may not be used without permission.