

CLPOSTSORT4000

POST ANALYTICS CLINICAL LABORATORY AUTOMATION

Output 900 to 4000 PPH




Powered by the precision of SuperTrak™ CONVEYANCE, the CLPostSort4000 automates post analytic processing of specimens following analytics in high throughput clinical labs from 900 to 4000 parts per hour.

Our modular approach allows us to deliver customized processing lines using standardized components, allowing you to easily adapt to evolving workflow requirements and scale your laboratory operations without disruption.

- **REALTIME MONITORING:** Illuminate™ Manufacturing Intelligence works with existing LIMS. Efficient recall processes and reporting with audit trails ensure transparency and accountability
- **OPERATIONAL EFFICIENCY:** Highest output per square foot; configurable systems designed to fit seamlessly into existing lab layouts
- **WORKFORCE OPTIMIZATION:** Reduce reliance on manual labor and sustain performance during staffing challenges by freeing up skilled personnel for higher-value tasks
- **PERFORMANCE & RELIABILITY:** Minimize human error and touchpoints with consistent, high-quality specimen handling

Whether you're dealing with staffing shortages or scaling to meet demand, ATS helps you increase lab productivity, reduce errors, and accelerate reliable results.


**EASILY AUTOMATE POST ANALYTIC PROCESSES FOR GREATER
EFFICIENCY, COMPLIANCE, RETENTION AND STORAGE**



CLPOSTSORT4000

TURNKEY LAB AUTOMATION SOLUTIONS STREAMLINE POST-ANALYTIC WORKFLOWS—EMPOWERING LABS TO DO MORE WITH LESS


Modular design allows you to arrange layout to accommodate floor space, lab requirements, throughput and turnaround time.



- 1 Input sorters x4 2 Sorter handoff x4 3 Main SuperTrak™
- 4 Sealer x10 5 Sorter handoff x4 6 Output sorters x4



Recapping



Sorting

POST-ANALYTIC PROCESSING

Following analysis, our systems ensure secure handling, storage, and disposal of specimens—supporting compliance and operational efficiency.

TAILORED OUTPUT SOLUTIONS



RECAPPING / RESEALING
Plug style capping or foil sealing



ARCHIVING
Output decks for sortation of archive targets into racks for organized storage



SAMPLE RETENTION & STORAGE
Configurable retention periods with organized archival systems



DISPOSAL AUTOMATION
Hands-free hazmat removal of expired or processed specimen



WORKFORCE OPTIMIZATION
Reduce reliance on manual labor; minimize human error with consistent, high-quality specimen handling



SCALABLE GROWTH
Expand capacity and capabilities without costly facility upgrades

READY TO FUTURE-PROOF YOUR LAB INFRASTRUCTURE?

- Modular, configurable systems designed to fit seamlessly into existing lab layouts and grow with your lab's evolving needs
- Maintain accuracy and control across every step of the analytic process

Whether upgrading legacy systems or launching a new lab, ATS Life Sciences Systems provides the proven automation solutions to keep your operations running smoothly, efficiently, and reliably.

Contact us today to schedule a demo or speak with a solutions expert.



[ATSLifeSciences.com](https://www.ATSLifeSciences.com) | lifesciences@atsautomation.com

ATS and the ATS LOGO are trademarks of ATS Corporation and may not be used without permission.